

## Soft Transfer Belt



- The Soft Transfer Belt is a popular belt used for transferring patients
- Two Carers normally work together when using this product. If necessary, one Carer can work alone, which depends on the patient and Carer's ability. If used in a care facility it may also depend on the facility's manual handling policy
- The belt is padded for comfort. The inside is a soft velour material. The outside is a nylon material
- There is an adjustable buckle fastening which adjusts easily to fit a variety of people
- The adjustable buckle allows you to keep it in the middle of the person when fastened
- It is suitable when the patient has a pendulous bust or if their body shape requires a narrower belt fastening
- The soft padded belt is slip resistant
- There are four comfortable, padded handles

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Changing body position

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Colour
TSE672240	Soft Transfer Belt	800 mm - 1600 mm	155 mm	Blue