

## Caring Cutlery



- A range of stainless steel utensils moulded into contoured, ivory-coloured plastic handles with built-up ends. These are of value to those with a weak grip and restricted wrist or finger movement
- The knife and fork have a shaped indent on the top for the index finger, to help with directional control
- The knife also has a serrated blade to assist with cutting. In addition to the standard spoon, there is a right handed angled spoon with an angled handle to minimise wrist movement
- Hygienically sealed, easy to clean and dishwasher safe
- Handle Length: 127 mm
- Weight: 50 - 65 g

## Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Eating and drinking

**Category:** Low risk

Click [here](#) for more information

### NDIS

**Item Code:**

05\_150900319\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Unit Weight
DLK276300	Knife	50 g
DLK276310	Fork	50 g
DLK276320	Spoon	60 g

Code	Display Name	Unit Weight
DLK276330	Caring Angled Spoon Right Hand	50 g