

## Solaris - ReadyWrap - Gauntlet



- Created to replace daytime compression bandaging or serve as an alternative/supplement to compression garments
- Inelastic materials produce high working pressure and low resting pressure
- Breathable and wickable materials ensure maximum comfort
- Provides coverage from the wrist to the MCP and overlaps the ReadyWrap arm garment to reduce gapping
- Left / Right specific design contours more naturally to the hand
- Pocket design allows for Swell Spots or foam pads to provide additional support

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Total Length	Colour	Wrist Circumference
FBS303659	Solaris Ready Wrap - Gauntlet - LEFT Large - Beige	145 mm	Beige	175 - 230 mm
FBS303660	Solaris Ready Wrap - Gauntlet - LEFT Large - Black	145 mm	Black	175 - 230 mm
FBS303661	Solaris Ready Wrap - Gauntlet - LEFT Medium - Beige	135 mm	Beige	155 - 200 mm
FBS303662	Solaris Ready Wrap - Gauntlet - LEFT Medium - Black	135 mm	Black	155 - 200 mm
FBS303663	Solaris Ready Wrap - Gauntlet - LEFT Small - Beige	125 mm	Beige	135 - 170 mm
FBS303667	Solaris Ready Wrap - Gauntlet - RIGHT Large - Beige	145 mm	Beige	175 - 230 mm
FBS303668	Solaris Ready Wrap - Gauntlet - RIGHT Large - Black	145 mm	Black	175 - 230 mm

Code	Display Name	Total Length	Colour	Wrist Circumference
FBS303669	Solaris Ready Wrap - Gauntlet - RIGHT Medium - Beige	135 mm	Beige	155 - 200 mm
FBS303670	Solaris Ready Wrap - Gauntlet - RIGHT Medium - Black	135 mm	Black	155 - 200 mm
FBS303672	Solaris Ready Wrap - Gauntlet - RIGHT Small - Black	125 mm	Black	135 - 170 mm