

## Aspire ComfiMotion Relax Bed



The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Box-on-Box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **USB Charge:** Connect and charge your devices conveniently.
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- **Leg Raise:** Promotes circulation, relieving tired, swollen and/or painful limbs.
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button.
- Total height adjustable in increments of 320, 370, 380 & 430 mm
- Frame - 10 Years
- Electrics - 2 Years
- Handset - 1 Year



Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046900	Aspire ComfiMotion Relax Bed - Portobello Grey - Long Single (Headboard Sep)	2030 mm	920 mm	320 - 430 mm	180 kg
BEB046910	Aspire ComfiMotion Relax Bed - Portobello Grey - King Single (Headboard Sep)	2030 mm	1070 mm	320 - 430 mm	180 kg
BEB046920	Aspire ComfiMotion Relax Bed - Portobello Grey - Long Double (Headboard Sep)	2030 mm	1370 mm	320 - 430 mm	360 kg
BEB046930	Aspire ComfiMotion Relax Bed - Portobello Grey - Long Queen (Headboard Sep)	2030 mm	1520 mm	320 - 430 mm	360 kg