


Monaco Chair



- The Monaco is ideal for users with reduced mobility or sitting balance who require comfort and support without compromising independence in transfers and function
- Its unique angle adjustable seat creates an optimum and safe sitting position when adjusted posteriorly. It reduces sliding, discomfort and shear as it helps to maintain positioning of the pelvis. It can also be adjusted to assist with stand transfers and independent ambulation
- Together with adjustable back angle recline, waterfall back, optional lateral supports, and pressure management materials, a high level of positioning and pressure redistribution is achievable without compromising the function or independence of the user
- The Monaco accommodates those with tight hamstrings as the negative angle clearance underneath the seat is designed to facilitate safe retraction of the feet behind the knees, ensuring the user's positioning is not dangerously disrupted
- The optional variable angle leg rest can be used to elevate the lower limbs and assist in the management of oedema. The inclusion of swivel castors, push handle and sliding footplate facilitate easy mobility and access to other environments within the home or care facility

Warranty Information

- 7 years frame
- 1 year upholstery

 Contract approved

ndis

Item Code:

05_180921181_0103_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Code	Display Name	Seat Depth	Seat Width	Total Height	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight
CHP197827	Chair Pressure Relief - Monaco Chair - Large - Manual - 160kg	350 - 550 mm	380 - 600 mm	1205 mm	800 mm	350 - 500 mm	160 kg	57 kg