


Aspire ComfiMotion Plush Pillow



- Plush Fibrefill sleep surface and Memory Foam core provides ergonomic support and comfort
- **Luxurious Fibrefill Outer:** Clouds of soft fibrefill provide a similar sleep experience to a traditional plush pillow
- **Ventilated Memory Foam Core:** Provides side and back sleepers with enhanced postural support, comfort and breathability
- **Breathable Easy-Zip Cover:** Encased in a removable cover for fuss-free care
- Familiar sleep surface
- Improved head and neck support
- Accomodates side and back sleepers
- Luxurious Plush Layer
- Immersive Memory Foam
- Medium Density
- Removable Cover

 Contract approved

Warranty Information

- 1 Year

| Code | Display Name | RRP | Box Quantity | Total Length | Total Width | Total Height |
|-----------|---------------------------------|-----|--------------|--------------|-------------|--------------|
| BEP050110 | Aspire ComfiMotion Plush Pillow | 128 | 1 | 660 mm | 480 mm | 110 mm |