

## Aspire ComfiMotion Plush Pillow



- Plush Fibrefill sleep surface and Memory Foam core provides ergonomic support and comfort
- **Luxurious Fibrefill Outer:** Clouds of soft fibrefill provide a similar sleep experience to a traditional plush pillow
- **Ventilated Memory Foam Core:** Provides side and back sleepers with enhanced postural support, comfort and breathability
- **Breathable Easy-Zip Cover:** Encased in a removable cover for fuss-free care
- Familiar sleep surface
- Improved head and neck support
- Accomodates side and back sleepers
- Luxurious Plush Layer
- Immersive Memory Foam
- Medium Density
- Removable Cover

✔ Contract approved

### Support @ Home

**Function:** Body stabilisation

**Category:** Under advice

Click [here](#) for more information

### Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height
BEP050110	Aspire ComfiMotion Plush Pillow	660 mm	480 mm	110 mm