


## Promedicare Inserto Tanto



- Foam Cushion designed for those who need a medium/high postural support
- Ideal for individuals with mild leg length discrepancies, wind-swept, pelvic obliquity, pelvic rotations, posterior and anterior pelvis tilt
- 13 modular and configurable inserts as standard
- Flexible open celled foam with low spring back action
- Comes with washable cover

 Contract approved

**ndis**

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:**05\_429\_0103\_1\_2

Click [here](#) for more information

### Warranty Information

- 2 Years Cushion
- 6 Months Cover



Code	Display Name	Safe Working Load	Pressure Care Rating
PTC608170	Promedicare Inserto Tanto - 18in x 20in	Unlimited	Medium Risk
PTC608172	Promedicare Inserto Tanto - 19in x 18in	Unlimited	Medium Risk