

## Aspire Vogue Indoor Walker



Ideal for use within the home, the Indoor Walker is perfect for portable dining and room to room mobility support.

- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Now supplied with bonus cane holder (ref. image)
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **No lifting & spillage** – meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable.
- **Sitting down for a meal** – when the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal.
- **Unique folding design** – the slimline design reduces the footprint when folded but allows the folded walker to stand upright on its own.
- **Brakes** – the ergonomic 'Soft Touch' lockable design allows the walker to be secured in place with minimal hand pressure.
- **Removable storage bag** – utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house.
- \*Unit weight includes accessories, bag and tray
- \*\*Frame weight is minimal functional weight

### Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Walking, manipulated by both arms

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_120606096\_0105\_1\_2

**Walking Aids -**

**Composite:**05\_416\_0105\_1\_2

Click [here](#) for more information

### Compliance

ISO 11199-2:2005

### Accessories

Click on the cart icon to display the range of optional accessories

Code	Display Name	Total Height	Safe Working Load	Unit Weight	Frame Weight	Colour
WAF705100CH	Aspire Vogue Indoor Walker - Champagne/Black	860 - 970 mm	110 kg	5.1 kg	4.4 kg	Champagne
WAF705100RE	Aspire Vogue Indoor Walker - Red/Black	860 - 970 mm	110 kg	5.1 kg	4.4 kg	Red