

Liko SupportVest (Mod. 91)



- Comfortable sling which lifts behind the back and under the arms
- Easy to use and requires only minimal individual adjustment so suit any patient
- Soft padding redistributes pressure evenly
- Supportive belt holds the vest in place at all times
- Available in Polyester

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394120	Support Vest Sling - Padding - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394130	Support Vest Sling - Padding - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394140	Support Vest Sling - Padding - Polyester - LARGE	200 kg	Polyester (Loops)

Code	Display Name	Safe Working Load	Material
LSS394150	Support Vest Sling - Padding - Polyester - XLARGE	200 kg	Polyester (Loops)