


Bed Exit Alarm



- Bed Exit Alarm Kit includes Mat, adaptor, AC plug & monitor
- Bed exit alarm fits on a bed mattress about kidney level
- When a patient gets out of the bed, the nurse call alarm is triggered
- Durable mat with tapered edges
- Suits most nurse call systems

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_222704219_0123_1_2

Personal Alarm – Composite:

05_427_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

| Code | Display Name | Total Length | Total Width |
|-----------|--|--------------|-------------|
| BEA007512 | Bed Exit Alarm Kit - Suits Nurse Call System | 760 mm | 250 mm |