

## Bed Exit Alarm



- Bed Exit Alarm Kit includes Mat, adaptor, AC plug & monitor
- Bed exit alarm fits on a bed mattress about kidney level
- When a patient gets out of the bed, the nurse call alarm is triggered
- Durable mat with tapered edges
- Suits most nurse call systems

### Warranty Information

- 1 Year



**ndis**

**Item Code:**

05\_222704219\_0123\_1\_2

**Personal Alarm – Composite:**

05\_427\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
BEA007512	Bed Exit Alarm Kit - Suits Nurse Call System	760 mm	250 mm