

## Freedom Hybrid Transroller Seat Walker/Wheelchair - BRO199



- Multi- use seat walker, rollator, wheelchair
- Light touch brakes
- Comfortable sling back rest (conforms to the shape of the user)
- Height adjustable handles
- Woven nylon seat
- Removable shoulder bag
- Large castors suitable for indoor or outdoor use
- Folds for easy transportation
- Double-X Frame for additional strength
- Swing away/ removable legrests
- Folded width: 270 mm

 Contract approved

### Support @ Home

**Function:** Walking, manipulated by both arms

**Category:** Prescribed

Click [here](#) for more information

### Compliance

ISO 11199-2:2005



### Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height
WAF750340	Freedom Hybrid Transroller Seat Walker/Wheelchair - Silver	710 - 980 mm	640 mm	860 - 980 mm	136 kg	9.3 kg	Silver	460 mm	560 mm