

Freedom Hybrid Transroller Seat Walker/Wheelchair - BRO199





- Multi- use seat walker, rollator, wheelchair
- Light touch brakes
- Comfortable sling back rest (conforms to the shape of the user)
- Height adjustable handles
- Woven nylon seat
- · Removable shoulder bag
- Large castors suitable for indoor or outdoor use
- Folds for easy transportation
- Double-X Frame for additional strength
- Swing away/ removable legrests
- Folded width: 270 mm

Warranty Information

• 1 Year

\odot		
	Contract	approved

Compliance

ISO 11199-2:2005

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height
WAF750340	Freedom Hybrid Transroller Seat Walker/Wheelchair - Silver	710 - 980 mm	640 mm	860 - 980 mm	136 kg	9.3 kg	Silver	460 mm	560 mm