

GelBodies Knee Protector



- Combined with pressure relieving gel, which is found in high-end wheelchair cushions, to offer additional padding to points on the body where pressure sores tend to form
- Improvement of circulation provides even distribution of the body weight and gives patient added overall comfort helping them to relax
- Breathable, stretchy cool-band material
- Sizing: Measure circumference above the knee. Choose the larger size if measurement is between sizes

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533325	GelBodies - Knee Protector - Sml	23-28cm
PTA533325M	GelBodies - Knee Protector - Med	28-33cm
PTA533325L	GelBodies - Knee Protector - Lge	33-41cm