

DONJOY AIRCAST AC-05 PNEUMATIC TENNIS ELBOW SUPPORT WITH AIR PAD



- Focused compression for less constriction.
- Breathable material enhances comfort and wearability.
- Simple, effective design makes it easy to apply.
- One size fits most for reduced inventory requirements. Application:
- Tennis elbow
- Golfer's Elbow

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code

Display Name

FBS306070 DONJOY AIRCAST AC-05 PNEUMATIC TENNIS ELBOW SUPPORT WITH AIR PAD - UNIVERSAL