

## Meiko Pan Rack



- Drying and storage rack for care utensils
- Stainless steel
- Drip tray is sold separately, size is compatible with storage rack

 Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information



Code	Display Name	Total Length	Total Width	Total Height	Capacity
UPA400300	Meiko Pan Rack - TS100 - 12 Pans	350 mm	1000 mm	650 mm	12 Pans
UPA400310	Meiko Pan Rack - TS50 - 3 Pans	350 mm	500 mm	650 mm	3 Pans
UPA400320	Meiko Pan Rack - TS70 - 6 Pans	350 mm	700 mm	650 mm	6 Pans
UPA400330	Meiko Pan Rack - TS90 - 9 Pans	350 mm	900 mm	650 mm	9 Pans
UPA400340	Meiko Pan Rack - 4 Pans - For vertical storage of bedpans (not mounted)	272 mm	444 mm	287 mm	4 Pans
UPA400350	Meiko Pan Rack - 8 Pans - For vertical storage of bedpans (not mounted)	272 mm	884 mm	287 mm	8 Pans