

Wall Cabinet with Combination Lock



- A wall mounted lockable chemical storage cabinet with the capacity for 2 x 5 litre bottles of either EcoWash+ or EcoCleanse Ultimate chemicals
- The cabinet is supplied with adaptor and flexible conduit to support the chemical supply tube linked directly to the machine
- Designed to hold up to 8 x 1 litre Click & Fit chemical bottles for the exclusive Pulpmatic Eco+ chemical dosing system
- Comes with combination lock

 Contract approved

Support @ Home

Function: Collecting urine and faeces

Category: Prescribed

Click [here](#) for more information

Compliance

ISO 15883/AS 4187:2019



Code	Display Name
UPA300300	Wall Cabinet - Combination Lock