

IBICI 'Repomen' Mens Compression Stockings



IBICI Repomen (for Men)

- Compression: 15 - 20mmHg
- Stylish appearance of dress socks
- Colours: Black, Beige, White

IBICI general overview:

Designed and manufactured in Italy, IBICI's compression legwear range is ideal for mild vascular conditions, or as a preventative measure. All IBICI garments have a closed toe (full foot). The range consists of:

- Segreta 140 (for Women)
- Segreta Plus (for Women)
- Repomen (for Men)
- Skywalk (Unisex)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641888	Ibici Repomen - Knee - 15-20mmHg - C/Toe - X Sml - Bge	Beige
PTS641888BEL	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Lge - Bge	Beige
PTS641888BEM	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Med - Bge	Beige
PTS641888BES	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Sml - Bge	Beige
PTS641888BEXL	Ibici Repomen - Knee - 15-20mmHg - C/Toe - XL - Bge	Beige

Code	Display Name	Colour
PTS641888BLL	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Lge - Blk	Black
PTS641888BLM	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Med - Blk	Black
PTS641888BLS	Ibici Repomen - Knee - 15-20mmHg - C/Toe - Sml - Blk	Black
PTS641888BLXL	Ibici Repomen - Knee - 15-20mmHg - C/Toe - XL - Blk	Black
PTS641888BLXS	Ibici Repomen - Knee - 15-20mmHg - C/Toe - X Sml - Blk	Black