

Aspire High Back Classic Day Chair



- A high back orthopaedic chair that offers adjustability, comfort and support
- Contoured Backrest – High back design provides increased postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Sturdy Design - Generous 160 kg weight capacity, yet lightweight for easy manoeuvrability
- Powder Coated Steel Frame – Easy to clean, strong and durable
- Non-Slip Tips – Provide stability whilst sitting or transferring
- Available in either Slate or Champagne Vinyl

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

NDIS

Item Code:

05_180915811_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 4688.2:2000; AS/NZS

4688.3:2000

Code	Display Name	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Backrest Height
CHP208050	Aspire High Back Classic Day Chair - Champagne Vinyl	160 kg	Steel, foam, vinyl	Champagne Vinyl	490 mm	490 mm	425 - 610 mm	720 mm

Code	Display Name	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Backrest Height
CHP208055	Aspire High Back Classic Day Chair - Slate Vinyl	160 kg	Steel, foam, vinyl	Slate Vinyl	490 mm	490 mm	425 - 610 mm	720 mm