

Aspire High Back Classic Day Chair





- A high back orthopaedic chair that offers adjustability, comfort and support
- Contoured Backrest High back design provides increased postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Easy clean vinyl surfaces
- Height Adjustable Easy push button adjustment for comfortable positioning and changing needs
- Sturdy Design Generous 160 kg weight capacity, yet lightweight for easy manoeuvrability
- Powder Coated Steel Frame Easy to clean, strong and durable
- Non-Slip Tips Provide stability whilst sitting or transferring
- Available in either Slate or Champagne Vinyl

Warranty Information

• 1 Year

- Contracted DVA product
- Ocontract approved

ndis

Item Code:

05_188800051_0103_1_2 Click **here** for more information

Compliance

AS/NZS 4688.2:2000; AS/NZS 4688.3:2000

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP208050	Aspire High Back Classic Day Chair - Champagne Vinyl	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Champagne Vinyl
CHP208055	Aspire High Back Classic Day Chair - Slate Vinyl	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Slate Vinyl