

Ottobock Backster Backrest



- Easily demountable one-hand release-back system
- Individually adaptable, five-piece padding with velcro connection to suit user requirements
- Stable back shell made of aluminium
- Can be mounted on nearly all commercially available wheelchairs
- Available in 3 key back-width sizes



ndis

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_417_0105_1_2

05_419_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years Backrest
- 6 Months Cover



Code	Display Name	Seat Width	Backrest Width	Backrest Height
PTC538000	Ottobock Backster Backrest 38 cm	380 - 420 mm	380 mm	450 mm
PTC538020	Ottobock Backster Backrest 42 cm	420 - 460 mm	420 mm	500 mm
PTC538010	Ottobock Backster Backrest 46 cm	460 - 500 mm	460 mm	550 mm