

Aidacare Emergency Cart



- The Aidacare Emergency Crash Cart has sleek lines, and strong construction
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Large work area with slide out work surface are standard
- Easy roll 125 mm castors
- Smooth manoeuvrability will make it the choice in emergency departments
- The simple breakaway lock offers easy access, combined with essential security
- Carts come standard with all essential accessories, to enhance the user experience.
- **Standard Inclusions:** Defibrillator Shelf, Cardiac Chest Board, Oxygen Bottle Holder, Extended Sliding Surface, Waste Bin.
- Customisation to suit a range of applications are also available upon request
- Comprehensive range of standard accessories available
- Standard drawer configuration detailed below
- **Small Emergency Crash Cart:**
 - 3 small 3' drawers (7 cm)
 - 1 medium 6' drawer (16 cm)
 - 1 large 9' drawer (25 cm)
- **Medium Emergency Crash Cart:**
 - 2 small 3' drawers (7 cm)
 - 2 medium 6' drawers (16 cm)
 - 1 large 9' drawer (25 cm)
- **Large Emergency Crash Cart:**
 - 3 small 3' drawers (7 cm)
 - 2 medium 6' drawers (16 cm)



Accessories

Click on the cart icon to display the range of optional accessories

- 1 large 9' drawer (25 cm)

Warranty Information

- Structural: 5 Years Warranty
- Accessories: 1 Year Warranty

| Code | Display Name | Total Length | Total Width | Total Height | Colour |
|-----------|----------------------------------|--------------|-------------|--------------|--------|
| WAT909400 | Aidacare Emergency Cart - Small | 790 mm | 610 mm | 900 mm | Red |
| WAT909410 | Aidacare Emergency Cart - Medium | 790 mm | 610 mm | 1000 mm | Red |
| WAT909420 | Aidacare Emergency Cart - Large | 790 mm | 610 mm | 1100 mm | Red |