

Aspire Deluxe Hygiene Access Sling - Polyester



- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding, quilted leg straps and a reinforced wide-cut back for comfort. The chest strap features velcro and buckle fastening for easy positioning and added safety
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- User Safety - smart feature such as velcro/buckle combinations for simplified positioning and greater user safety
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 2 or 3

Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_123603163_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	RRP	Box Quantity	Safe Working Load	Material
LSS390760	Deluxe Hygiene Access Sling - Polyester - Small	392	1	200 kg	Polyester
LSS390765	Deluxe Hygiene Access Sling - Polyester - Medium	392	1	200 kg	Polyester
LSS390770	Deluxe Hygiene Access Sling - Polyester - Large	392	1	300 kg	Polyester
LSS390775	Deluxe Hygiene Access Sling - Polyester - Extra Large	392	1	300 kg	Polyester
LSS390776	Deluxe Hygiene Access Sling - Polyester - XXL	420	1	300 kg	Polyester (Loops)
LSS390777	Deluxe Hygiene Access Sling - Polyester - XXXL	433	1	300 kg	Polyester (Loops)