

Dysphagia Cup



- Designed to assist people who have trouble swallowing, drinking while tucking the chin toward the chest gives normal swallowing mechanisms time to work
- The design of this cup helps prevent liquids from escaping at the lips and directs the liquid to the centre of the mouth
- Weight: 230 g

Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Capacity
DLK275470	Dysphagia Cup - Green	200 ml