


## Oxford Journey Folding Standing Lifter



- Mobile compact patient lifter for transfer to a seated, lying or standing position
- Compact and lightweight allowing carers to manoeuvre in tight spaces
- Folds away for easy storage or transportation
- Unique adjustable 'cow-horn' shape supports a wide range of patient heights and sizes
- Adjustable knee pad
- Epoxy coated frame
- Electric hand control operation
- Emergency operation
- Easily manoeuvrable with 125 mm rear locking castors and 75 mm front twin swivel castors
- Electric leg gating
- Supplied with battery and charger
- Item Weight: 40.2 kg

 Contracted DVA product

 Contract approved

**ndis**

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS377305	Journey Folding Standing Lifter	1000 mm	685 mm	1130 - 1600 mm	720 - 1530 mm	590 - 875 mm	True	155 kg	Aluminium