

Oxford Journey Folding Standing Lifter



- Mobile compact patient lifter for transfer to a seated, lying or standing position
- Compact and lightweight allowing carers to manoeuvre in tight spaces
- Folds away for easy storage or transportation
- Unique adjustable 'cow-horn' shape supports a wide range of patient heights and sizes
- Adjustable knee pad
- Epoxy coated frame
- Electric hand control operation
- Emergency operation
- Easily manoeuvrable with 125 mm rear locking castors and 75 mm front twin swivel castors
- Electric leg gating
- Supplied with battery and charger
- Item Weight: 40.2 kg

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS377305	Journey Folding Standing	1000 mm	685 mm	1130 - 1600	155 kg	Aluminium	40.2 kg	720 - 1530	590 - 875 mm	True

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
	Lifter			mm				mm		