

Spex Flex Cushion



- Excellent choice for individuals who have symmetrical or asymmetrical limited hip flexion.
- Strategic positioning base, pads and wedges offer complete customisations on-the-spot for a perfect fit for the user.
- Bevelled front edge allows for knee contractures
- Positioning base primarily accommodates lack of hip flexion up to 20° but due to its innovative positioning base it can support ruling out other associated deformities as pelvic obliquity, windsweeping, hip abduction, hip internal rotation with tight hamstrings and even lack of hip extension

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name
------	--------------

PTC607284	SPEX Flex System - 16x16in
-----------	----------------------------