

Jazzy 1450



The Jazzy® 1450 is a front-wheel drive bariatric power base which has a weight capacity of 272 kg (600 lbs.) and durable construction to make it an exceptionally rugged power chair. Its unique design delivers high-performance operation indoors and out.

Warranty Information

- Frame: Lifetime limited
- Electronics and Drive Motors: 13 Months



 Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122306125_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat Depth	Seat Width	Maximum Speed
PWS654269	Jazzy 1450	970 mm	742 mm	273 kg	812 mm	609 mm	7.2 km/h