

## Aspire Homecare Over Toilet Aid



- Adjustable height raised toilet seat allows for safe ingress and egress when toileting
- Durable high grip feet for stability and safety
- Adjustable legs to suit a wide range of users
- Raised armrests assist the user in transferring on and off
- Width between the arms: 475 mm

 Contract approved

### Support @ Home

**Function:** Toileting

**Category:** Prescribed

Click [here](#) for more information

### Warranty Information

- 1 Year



Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTT145025	Aspire Homecare Over Toilet Aid	520 mm	595 - 750 mm	100 kg	410 mm	340 mm	405 - 560 mm	405 - 560 mm