

## Solaris - ReadyWrap - Knee



- ReadyWrap garments are designed to aid people living with venous disease or insufficiencies
- Created to replace daytime compression bandaging or serve as an alternative or supplement to elastic compression stockings
- Silver microfibre under-garment liner is available for added comfort
- Contraindications: arterial insufficiency or degeneration, deep vein thrombosis, untreated congestive heart failure, localised or systemic infection or cancer, absent or severely impaired sensation
- The ReadyWrap knee unit provides support for those with localised patella swelling or acts as a bridge connecting the thigh and calf units to provide full leg coverage
- The simple 3-strap design makes that unit easy to don and the 50% overlap mimics short-stretch bandaging

☑ Contracted DVA product

☑ Contract approved

Code	Display Name	Colour
FBS303629	Solaris ReadyWrap - Knee - Large - Beige	Beige
FBS303630	Solaris ReadyWrap - Knee - Large - Black	Black
FBS303631	Solaris ReadyWrap - Knee - Medium - Beige	Beige
FBS303632	Solaris ReadyWrap - Knee - Medium - Black	Black
FBS303633	Solaris ReadyWrap - Knee - Small - Beige	Beige
FBS303634	Solaris ReadyWrap - Knee - Small - Black	Black

<b>Code</b>	<b>Display Name</b>	<b>Colour</b>
FBS303635	Solaris ReadyWrap - Knee - Xlarge - Beige	Beige
FBS303637	Solaris ReadyWrap - Knee - XXlarge - Beige	Beige
FBS303638	Solaris ReadyWrap - Knee - XXlarge - Black	Black