

Bedside Commode Maxi Fixed Height



- Large bedside chair can be used for toileting
- Extra seat width appropriate for maximum user weight of up to 210 kg
- Fold back padded seat for easy access to commode bowl
- Contoured backrest for added comfort and support
- Comes with lift up toilet seat flap, fixed bowl carrier and bowl with removable lid and carry handle
- Fitted with non-slip rubber tips
- Reinforced mild steel epoxy coated frame is zinc treated for corrosion resistance
- Champagne colour vinyl



Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
BEC046040	Maxi Commode - Champagne - 600mm Wide	195 kg	15 kg	Champagne Vinyl	490 mm	600 mm	500 mm
BEC046050	Maxi Commode - Champagne - 650mm Wide	210 kg	15 kg	Champagne Vinyl	490 mm	650 mm	500 mm