

## Oppo 2031 Knee Open Patella Support



- May be used for stabilisation and support of the knee following strains, injury or surgery
- Elastic brace with motion hinges provides medial and lateral support for the unstable knee
- Provides firm support for the maintenance of joint structures with gentle compression over the kneecap
- For arthritic knee conditions, minor sprains and strains
- Two hook and loop straps insure comfortable fit

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Lower limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS301530M	Oppo - Knee Hinged - Pull-on Elastic - O/Patella - Med	343 - 381 mm
FBS301530L	Oppo - Knee Hinged - Pull-on Elastic - O/Patella - Lge	381 - 419 mm
FBS301530	Oppo - Knee Hinged - Pull-on Elastic - O/Patella - Sml	305 - 343 mm
FBS301530XL	Oppo - Knee Hinged - Pull-on Elastic - O/Patella - XL	419 - 457 mm